



U-ramps for manual wheelchairs Telescopic and Fixed PERFO-LIGHT ramps

Our ramps are light-weight, easy to handle and just the right width.

perfo
L I G H T

Nylon glides ensure the telescopic ramps extend smoothly and virtually silently. A locking device prevents them sliding part.



Everyday, wheelchair users encounter many varied step heights, needing different ramp lengths. Telescopic ramps enable you to choose the suitable length, taking into account available space and/or gradient required.

Ideal for wheelchair users.



V-ramps Foldable PERFO-LIGHT ramps

V-ramps are particularly suitable for manual wheelchairs, electric 3-wheel mobility scooters and can even be used for many electric 4-wheel wheelchair models.

perfo.
L I G H T

The width of the driving area is the same throughout the whole ramp and the inside of the kerbs is completely smooth.

Our foldable ramps each have two handles to keep the ramps together when you carry them.




These ramps compliment the telescopic ramp range. Some users find it easier to handle a folding ramp.

The upper ends of the side kerbs have a long bevel, ensuring that even a low clearance wheelchair is able to drive into e.g. a car, even though the angle may be steep.



FOR A 3-WHEEL MOBILITY SCOOTER YOU SIMPLY PLACE 3 RAMPS BESIDE EACH OTHER.



ONE ADVANTAGE WITH THE FOLDING RAMPS IS THAT THEY CAN BE USED LIKE A BRIDGE OVER A THRESHOLD OR DOORSTEP.

E-ramps

The E-ramp is a universal ramp, perfect for all wheelchairs

perfo.
L I G H T

The E-ramp is perfect for all wheelchairs, including 4 wheel electric wheelchairs which often require a wider driving width.


FOR WHEELCHAIRS THAT REQUIRE A WIDER DRIVING WIDTH



This ramp is highly recommended for all public buildings, including shops, restaurants and hotels.

See page 6 for technical information on all lengths and maximum heights for every ramp.

4 wheel electric wheelchairs often require a wider driving width






portable ramps

*PERFO-light offers
a wide range of
practical solutions*

Our PERFO-LIGHT Ramps offer a wide range of practical solutions to many challenges – and not just for wheelchair users.






They are also suitable for trolleys, motorbikes, for use at home or away, visiting friends and relatives.

Easy to use and easy to carry, the ramps come with handles that fold flat against the sides. Weight is reduced by perforating the ramps, providing an excellent anti-slip driving surface. The ramps can be set-up everywhere, without need for any tools.

PERFO LIGHT Ramps offer everyday solutions, safely and simply, providing greater freedom, ease and peace of mind.

General information

Art no.	Product name	Length min	Length max	internal width	Profile	Weight /each	Max load /pair
Telescopic ramps							
U 2	Perfolight telescopic	116 cm	202 cm	19,5 cm		5,6 kg	300 kg
U 2/3	Perfolight telescopic	78 cm	202 cm	18,0 cm		5,8 kg	300 kg
U 3	Perfolight telescopic	116 cm	288 cm	18,0 cm		8,2 kg	250 kg
Foldable and fixed ramps							
V 5	Perfolight pavement ramp	55 cm	55 cm	21,0 cm		1,7 kg	400 kg
V 10	Perfolight fixed	116 cm	116 cm	21,0 cm		3,1 kg	400 kg
V 11	Perfolight foldable	55 cm	110 cm	21,0 cm		3,6 kg	300 kg
V 15	Perfolight foldable	75 cm	150 cm	21,0 cm		4,5 kg	300 kg
V 20	Perfolight foldable	100 cm	200 cm	21,0 cm		5,9 kg	250 kg
Power wheelchair ramps							
E 1	Perfolight fixed	106 cm	106 cm	26,5 cm		3,7 kg	400 kg
E 2	Perfolight foldable	106 cm	212 cm	26,5 cm		7,4 kg	300 kg
E 2/3	Perfolight foldable, telescopic	72 cm	212 cm	25,0 cm		7,8 kg	300 kg
E 3	Perfolight foldable, telescopic	106 cm	300 cm	25,0 cm		11,0 kg	250 kg

Rampguide

Height:	10 cm	20 cm	30 cm	40 cm	60 cm
Ramp	V 5	V 10 V 11	V 15	U 2 U 2/3 V 20	U 3 E 3
		E 1		E 2	E 2/3

PLEASE NOTE: THIS IS ONLY A RECOMMENDATION.

Contact Information

Address:

Box 19007
SE-720 19 Västerås
Sweden
Phone: +46 21 661 00
Fax: +46 21 662 33

Visiting address:

Pistonggatan 4B
Västerås
Office hours: 7.30 am to 4.30 pm
info@kvisstbergprodukter.se
www.kvisstbergprodukter.se